

GDF-9 Polyclonal Antibody

Catalog No: #41706

Package Size: #41706-1 50ul #41706-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GDF-9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GDF-9 Polyclonal Antibody detects endogenous levels of GDF-9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human GDF-9.
Target Name	GDF-9
Other Names	GDF9; Growth/differentiation factor 9; GDF-9
Accession No.	Swiss-Prot: O60383NCBI Gene ID: 2661
Uniprot	O60383
GeneID	2661;
SDS-PAGE MW	51kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of L929 cells using GDF-9 Polyclonal Antibody

Note: This product is for in vitro research use only