

ARP Polyclonal Antibody

Catalog No: #41711

Package Size: #41711-1 50ul #41711-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ARP Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	ARP Polyclonal Antibody detects endogenous levels of ARP protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ARP.
Target Name	ARP
Other Names	MANF; ARMET; ARP; Mesencephalic astrocyte-derived neurotrophic factor; Arginine-rich protein; Protein ARMET
Accession No.	Swiss-Prot: P55145NCBI Gene ID: 7873
Uniprot	P55145
GeneID	7873;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using ARP Polyclonal Antibody

Note: This product is for in vitro research use only