

## Neuregulin-1 Polyclonal Antibody

Catalog No: #41715

Package Size: #41715-1 50ul #41715-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Neuregulin-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Neuregulin-1 Polyclonal Antibody detects endogenous levels of Neuregulin-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Neuregulin-1.
Target Name	Neuregulin-1
Other Names	NRG1; GGF; HGL; HRGA; NDF; SMDF; Pro-neuregulin-1, membrane-bound isoform; Pro-NRG1
Accession No.	Swiss-Prot: Q02297NCBI Gene ID: 3084
Uniprot	Q02297
GeneID	3084;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of HepG2 22RV1 HELA 293 AD293 cells using Neuregulin-1 Polyclonal Antibody

---

Note: This product is for in vitro research use only