

HNF-3 α / β / γ Polyclonal Antibody

Catalog No: #41718

Package Size: #41718-1 50ul #41718-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HNF-3 α / β / γ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	HNF-3 α / β / γ Polyclonal Antibody detects endogenous levels of HNF-3 α / β / γ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HNF-3 α / β / γ around the non-acetylation site of K265.
Target Name	HNF-3 α / β / γ
Other Names	FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta; HNF-3B; Forkhead box protein A2;
Accession No.	Swiss-Prot: P55317NCBI Gene ID: 3169
Uniprot	P55317
GeneID	3169;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of A549 cells using HNF-3 α / β / γ Polyclonal Antibody

Note: This product is for in vitro research use only