

HNF-4 α/γ Polyclonal Antibody

Catalog No: #41720

Package Size: #41720-1 50ul #41720-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HNF-4 α/γ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	HNF-4 α/γ Polyclonal Antibody detects endogenous levels of HNF-4 α/γ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HNF-4 α/γ around the non-acetylation site of K127.
Target Name	HNF-4 α/γ
Other Names	HNF4A; HNF4; NR2A1; TCF14; Hepatocyte nuclear factor 4-alpha; HNF-4-alpha; Nuclear receptor subfamily 2 group A member 1; Transcription factor 14; TCF-14;
Accession No.	Swiss-Prot: P41235NCBI Gene ID: 3172
Uniprot	P41235
GeneID	3172;
SDS-PAGE MW	52kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of Jurkat K562 cells using HNF-4 α/γ Polyclonal Antibody

Note: This product is for in vitro research use only