

MyoD Polyclonal Antibody

Catalog No: #41730

Package Size: #41730-1 50ul #41730-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	MyoD Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	MyoD Polyclonal Antibody detects endogenous levels of MyoD protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human MyoD around the non-acetylation site of K99/K102.
Target Name	MyoD
Other Names	MYOD1; BHLHC1; MYF3; MYOD; Myoblast determination protein 1; Class C basic helix-loop-helix protein 1; bHLHc1; Myogenic factor 3; Myf-3
Accession No.	Swiss-Prot: P15172NCBI Gene ID: 4654
Uniprot	P15172
GeneID	4654;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using MyoD Polyclonal Antibody

Note: This product is for in vitro research use only