

Nrf2 Polyclonal Antibody

Catalog No: #41731

Package Size: #41731-1 50ul #41731-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Nrf2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Nrf2 Polyclonal Antibody detects endogenous levels of Nrf2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Nrf2 around the non-acetylation site of K599.
Target Name	Nrf2
Other Names	NFE2L2; NRF2; Nuclear factor erythroid 2-related factor 2; NF-E2-related factor 2; NFE2-related factor 2; HEBP1; Nuclear factor, erythroid derived 2, like 2
Accession No.	Swiss-Prot: Q16236NCBI Gene ID: 4780
Uniprot	Q16236
GeneID	4780;
SDS-PAGE MW	67kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using Nrf2 Polyclonal Antibody

Note: This product is for in vitro research use only