

## MMP-9 Polyclonal Antibody

Catalog No: #41747

Package Size: #41747-1 50ul #41747-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	MMP-9 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	MMP-9 Polyclonal Antibody detects endogenous levels of MMP-9 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MMP-9.
Target Name	MMP-9
Other Names	MMP9; CLG4B; Matrix metalloproteinase-9; MMP-9; 92 kDa gelatinase; 92 kDa type IV collagenase; Gelatinase B; GELB
Accession No.	Swiss-Prot: P14780NCBI Gene ID: 4318
Uniprot	P14780
GeneID	4318;
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of K562 cells using MMP-9 Polyclonal Antibody.

---

Note: This product is for in vitro research use only