

## IGFBP6 Polyclonal Antibody

Catalog No: #41750

Package Size: #41750-1 50ul #41750-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | IGFBP6 Polyclonal Antibody                                                                                            |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB ELISA                                                                                                              |
| Species Reactivity    | Hu                                                                                                                    |
| Specificity           | IGFBP6 Polyclonal Antibody detects endogenous levels of IGFBP6 protein.                                               |
| Immunogen Type        | peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from the Internal region of human IGFBP6.                                                 |
| Target Name           | IGFBP6                                                                                                                |
| Other Names           | IGFBP6; IBP6; Insulin-like growth factor-binding protein 6; IBP-6; IGF-binding protein 6; IGFBP-6                     |
| Accession No.         | Swiss-Prot: P24592NCBI Gene ID: 3489                                                                                  |
| Uniprot               | P24592                                                                                                                |
| GeneID                | 3489;                                                                                                                 |
| SDS-PAGE MW           | 25kd                                                                                                                  |
| Concentration         | 1mg/ml                                                                                                                |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of K562 cells using IGFBP6 Polyclonal Antibody.

---

Note: This product is for in vitro research use only