

SCYL1BP1 Polyclonal Antibody

Catalog No: #41756

Package Size: #41756-1 50ul #41756-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	SCYL1BP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	SCYL1BP1 Polyclonal Antibody detects endogenous levels of SCYL1BP1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human SCYL1BP1.
Target Name	SCYL1BP1
Other Names	GORAB; NTKLBP1; SCYL1BP1; RAB6-interacting golgin; N-terminal kinase-like-binding protein 1; NTKL-BP1; NTKL-binding protein 1; hNTKL-BP1; SCY1-like 1-binding protein 1; SCYL1-BP1; SCYL1-binding protein 1
Accession No.	Swiss-Prot: Q5T7V8NCBI Gene ID: 92344
Uniprot	Q5T7V8
GeneID	92344;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of AD293 cells using SCYL1BP1 Polyclonal Antibody.

Note: This product is for in vitro research use only