

Annexin VII Polyclonal Antibody

Catalog No: #41760

Package Size: #41760-1 50ul #41760-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Annexin VII Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	Annexin VII Polyclonal Antibody detects endogenous levels of Annexin VII protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Annexin VII.
Target Name	Annexin VII
Other Names	ANXA7; ANX7; SNX; OK/SW-cl.95; Annexin A7; Annexin VII; Annexin-7; Synexin
Accession No.	Swiss-Prot: P20073NCBI Gene ID: 310
Uniprot	P20073
GeneID	310;
SDS-PAGE MW	53kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of A549 cells using Annexin VII Polyclonal Antibody.

Note: This product is for in vitro research use only