

CGREF1 Polyclonal Antibody

Catalog No: #41763

Package Size: #41763-1 50ul #41763-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CGREF1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CGREF1 Polyclonal Antibody detects endogenous levels of CGREF1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human CGREF1.
Target Name	CGREF1
Other Names	CGREF1; CGR11; Cell growth regulator with EF hand domain protein 1; Cell growth regulatory gene 11 protein; Hydrophobestin
Accession No.	Swiss-Prot: Q99674NCBI Gene ID: 10669
Uniprot	Q99674
GeneID	10669;
SDS-PAGE MW	32kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using CGREF1 Polyclonal Antibody.

Note: This product is for in vitro research use only