

Connexin-40 Polyclonal Antibody

Catalog No: #41764

Package Size: #41764-1 50ul #41764-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Connexin-40 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Connexin-40 Polyclonal Antibody detects endogenous levels of Connexin-40 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Connexin-40.
Target Name	Connexin-40
Other Names	GJA5; Gap junction alpha-5 protein; Connexin-40; Cx40
Accession No.	Swiss-Prot: P36382NCBI Gene ID: 2702
Uniprot	P36382
GeneID	2702;
SDS-PAGE MW	39kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using Connexin-40 Polyclonal Antibody.

Note: This product is for in vitro research use only