

## Apaf-1 Polyclonal Antibody

Catalog No: #41765

Package Size: #41765-1 50ul #41765-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Apaf-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Apaf-1 Polyclonal Antibody detects endogenous levels of Apaf-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Apaf-1.
Target Name	Apaf-1
Other Names	APAF1; KIAA0413; Apoptotic protease-activating factor 1; APAF-1
Accession No.	Swiss-Prot: O14727NCBI Gene ID: 317
Uniprot	O14727
GeneID	317;
SDS-PAGE MW	135kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of K562 cells using Apaf-1 Polyclonal Antibody.

---

Note: This product is for in vitro research use only