

IL-2R $\gamma$  Polyclonal Antibody

Catalog No: #41767

Package Size: #41767-1 50ul #41767-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	IL-2R $\gamma$ Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-2R $\gamma$ Polyclonal Antibody detects endogenous levels of IL-2R $\gamma$ protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-2R $\gamma$ .
Target Name	IL-2R $\gamma$
Other Names	IL2RG; Cytokine receptor common subunit gamma; Interleukin-2 receptor subunit gamma; IL-2 receptor subunit gamma; IL-2R subunit gamma; IL-2RG; gammaC; p64; CD132
Accession No.	Swiss-Prot: P31785NCBI Gene ID: 3561
Uniprot	P31785
GeneID	3561;
SDS-PAGE MW	40kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images

Western Blot analysis of K562 cells using IL-2R $\gamma$  Polyclonal Antibody.

---

Note: This product is for in vitro research use only