

17 β -HSD4 Polyclonal Antibody

Catalog No: #41772

Package Size: #41772-1 50ul #41772-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	17 β -HSD4 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	17 β -HSD4 Polyclonal Antibody detects endogenous levels of 17 β -HSD4 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human 17 β -HSD4.
Target Name	17 β -HSD4
Other Names	HSD17B4; EDH17B4; Peroxisomal multifunctional enzyme type 2; MFE-2; 17-beta-hydroxysteroid dehydrogenase 4; 17-beta-HSD 4; D-bifunctional protein; DBP; Multifunctional protein 2; MPF-2
Accession No.	Swiss-Prot: P51659NCBI Gene ID: 3295
Uniprot	P51659
GeneID	3295;
SDS-PAGE MW	80kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images

Western Blot analysis of HuvEc, K562 cells using 17 β -HSD4 Polyclonal Antibody.

Note: This product is for in vitro research use only