

CD42b Polyclonal Antibody

Catalog No: #41778

Package Size: #41778-1 50ul #41778-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	CD42b Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	CD42b Polyclonal Antibody detects endogenous levels of CD42b protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human CD42b.
Target Name	CD42b
Other Names	GP1BA; Platelet glycoprotein Ib alpha chain; GP-Ib alpha; GPIb-alpha; GPIbA; Glycoprotein Ibalpha; Antigen CD42b-alpha; CD42b
Accession No.	Swiss-Prot: P07359NCBI Gene ID:
Uniprot	P07359
GeneID	2811;
SDS-PAGE MW	69kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using CD42b Polyclonal Antibody.

Note: This product is for in vitro research use only