

Jagged1 Polyclonal Antibody

Catalog No: #41786

Package Size: #41786-1 50ul #41786-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Jagged1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IF IHC-p
Species Reactivity	Hu Ms Rt
Specificity	Jagged1 Polyclonal Antibody detects endogenous levels of Jagged1 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human Jagged1.
Target Name	Jagged1
Other Names	JAG1; JAGL1; Protein jagged-1; Jagged1; hJ1; CD339
Accession No.	Swiss-Prot: P78504NCBI Gene ID: 182
Uniprot	P78504
GeneID	182;
SDS-PAGE MW	140kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


ELISA: 1/20000.

Not yet tested in other applications.

IF: 1:50-200

IHC-p: 1/100-1/300

Images



Western Blot analysis of HeLa cells using Jagged1 Polyclonal Antibody.

Note: This product is for in vitro research use only