

## Human Regenerating gene 4 (REG-4) ELISA Kit

Catalog No: #EK7547



Package Size: #EK7547-1 48T #EK7547-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description

Product Name	Human Regenerating gene 4 (REG-4) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human (Homo sapiens)
Other Names	GISP; REG-IV; RELP; OTTHUMP00000014031 gastrointestinal secretory protein regenerating gene type IV
Accession No.	Q9BYZ8
Uniprot	Q9BYZ8
GeneID	83998;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

## Application Details

Detect Range:0.156-10 ng/mL

Sensitivity:0.057 ng/mL

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

**Detection Method:**SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate REG4 in samples. An antibody specific for REG4 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyREG4 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for REG4 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of REG4 bound in the initial step. The color development is stopped and the intensity of the color is measured.**Product Overview:**The REG4 gene encodes a deduced 158-amino acid protein with an N-terminal signal peptide of 22 amino acids and a predicted molecular mass of 18.2 kD. REG4 maintains conservation of the 6 cysteine residues present in other known members of the REG family. By Northern blot analysis of intestinal mucosa from 12 individuals with inflammatory bowel disease, Hartupee et al. (2001) identified a 1.1-kb REG4 transcript in all samples. Control specimens had low but detectable expression of REG4. REG4 expression was significantly upregulated by inflammation and tissue injury associated with active Crohn disease and ulcerative colitis. RNA dot-blot analysis showed that most of the tissues expressing REG4 are in the gastrointestinal tract, i.e., colon, small intestine, stomach, and pancreas.

---

Note: This product is for in vitro research use only