

EKLK Polyclonal Antibody

Catalog No: #41805

Package Size: #41805-1 50ul #41805-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	EKLK Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	EKLK Polyclonal Antibody detects endogenous levels of EKLK protein only when non-acetylation at K274.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human EKLK around the non-acetylation site of K274.
Target Name	EKLK
Other Names	KLF1; EKLK; Krueppel-like factor 1; Erythroid krueppel-like transcription factor; EKLK
Accession No.	Swiss-Prot: Q13351NCBI Gene ID: 10661
Uniprot	Q13351
GeneID	10661;
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of LOVO cells using EKLK Polyclonal Antibody.

Note: This product is for in vitro research use only