

c-Myb Polyclonal Antibody

Catalog No: #41806

Package Size: #41806-1 50ul #41806-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	c-Myb Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	c-Myb Polyclonal Antibody detects endogenous levels of c-Myb protein only when non-acetylation at K480.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human c-Myb around the non-acetylation site of K480.
Target Name	c-Myb
Other Names	MYB; Transcriptional activator Myb; Proto-oncogene c-Myb
Accession No.	Swiss-Prot: P10242NCBI Gene ID: 4602
Uniprot	P10242
GeneID	4602;
SDS-PAGE MW	70kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using c-Myb Polyclonal Antibody.

Note: This product is for in vitro research use only