

Tubulin α Polyclonal Antibody

Catalog No: #41811

Package Size: #41811-1 50ul #41811-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Tubulin α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Tubulin α Polyclonal Antibody detects endogenous levels of Tubulin α protein only when non-acetylation at K112.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Tubulin α around the non-acetylation site of K112.
Target Name	Tubulin α
Other Names	TUBA1A; TUBA3; Tubulin alpha-1A chain; Alpha-tubulin 3; Tubulin B-alpha-1; Tubulin alpha-3 chain; TUBA1B; Tubulin alpha-1B chain; Alpha-tubulin ubiquitous
Accession No.	Swiss-Prot: Q71U36/P68363/Q9BQE3NCBI Gene ID: 7846/10376/84790
Uniprot	Q71U36
GeneID	7846;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562, PC12 cells using Tubulin α Polyclonal Antibody.

Note: This product is for in vitro research use only