

GPR87/95 Polyclonal Antibody

Catalog No: #41816

Package Size: #41816-1 50ul #41816-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	GPR87/95 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	GPR87/95 Polyclonal Antibody detects endogenous levels of GPR87/95 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human GPR87/95.
Target Name	GPR87/95
Other Names	GPR87; GPR95; FKSG78; G-protein coupled receptor 87; G-protein coupled receptor 95
Accession No.	Swiss-Prot: Q9BY21NCBI Gene ID: 53836
Uniprot	Q9BY21
GeneID	53836;
SDS-PAGE MW	42kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using GPR87/95 Polyclonal Antibody.

Note: This product is for in vitro research use only