

ApoA-I Polyclonal Antibody

Catalog No: #41817

Package Size: #41817-1 50ul #41817-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ApoA-I Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms
Specificity	ApoA-I Polyclonal Antibody detects endogenous levels of ApoA-I protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ApoA-I.
Target Name	ApoA-I
Other Names	APOA1; Apolipoprotein A-I; Apo-AI; ApoA-I; Apolipoprotein A1
Accession No.	Swiss-Prot: P02647NCBI Gene ID: 335
Uniprot	P02647
GeneID	335;
SDS-PAGE MW	31kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using ApoA-I Polyclonal Antibody.

Note: This product is for in vitro research use only