

FGF-16 Polyclonal Antibody

Catalog No: #41818

Package Size: #41818-1 50ul #41818-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	FGF-16 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	FGF-16 Polyclonal Antibody detects endogenous levels of FGF-16 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human FGF-16.
Target Name	FGF-16
Other Names	FGF16; Fibroblast growth factor 16; FGF-16
Accession No.	Swiss-Prot: O43320NCBI Gene ID: 8823
Uniprot	O43320
GeneID	8823;
SDS-PAGE MW	24kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of K562 cells using FGF-16 Polyclonal Antibody.

Note: This product is for in vitro research use only