

ECA39 Polyclonal Antibody

Catalog No: #41820

Package Size: #41820-1 50ul #41820-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	ECA39 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	ECA39 Polyclonal Antibody detects endogenous levels of ECA39 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human ECA39.
Target Name	ECA39
Other Names	BCAT1; BCT1; ECA39; Branched-chain-amino-acid aminotransferase, cytosolic; BCAT(c); Protein ECA39
Accession No.	Swiss-Prot: P54687NCBI Gene ID: 586
Uniprot	P54687
GeneID	586;
SDS-PAGE MW	43kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of K562 cells using ECA39 Polyclonal Antibody.

Note: This product is for in vitro research use only