

MAPK Organizer 1 Polyclonal Antibody

Catalog No: #41821

Package Size: #41821-1 50ul #41821-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAPK Organizer 1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	MAPK Organizer 1 Polyclonal Antibody detects endogenous levels of MAPK Organizer 1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MAPK Organizer 1.
Target Name	MAPK Organizer 1
Other Names	WDR83; MORG1; WD repeat domain-containing protein 83; Mitogen-activated protein kinase organizer 1; MAPK organizer 1
Accession No.	Swiss-Prot: Q9BRX9NCBI Gene ID: 84292
Uniprot	Q9BRX9
GeneID	84292;
SDS-PAGE MW	34kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using MAPK Organizer 1 Polyclonal Antibody.

Note: This product is for in vitro research use only