

EP1 Polyclonal Antibody

Catalog No: #41823

Package Size: #41823-1 50ul #41823-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	EP1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu
Specificity	EP1 Polyclonal Antibody detects endogenous levels of EP1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human EP1.
Target Name	EP1
Other Names	PTGER1; Prostaglandin E2 receptor EP1 subtype; PGE receptor EP1 subtype; PGE2 receptor EP1 subtype; Prostanoid EP1 receptor
Accession No.	Swiss-Prot: P34995NCBI Gene ID: 5731
Uniprot	P34995
GeneID	5731;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of HepG2 cells using EP1 Polyclonal Antibody.

Note: This product is for in vitro research use only