

SDHB Polyclonal Antibody

Catalog No: #41824

Package Size: #41824-1 50ul #41824-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SDHB Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	SDHB Polyclonal Antibody detects endogenous levels of SDHB protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human SDHB.
Target Name	SDHB
Other Names	SDHB; SDH; SDH1; Succinate dehydrogenase [ubiquinone] iron-sulfur subunit, mitochondrial; Iron-sulfur subunit of complex II; Ip
Accession No.	Swiss-Prot: P21912NCBI Gene ID: 6390
Uniprot	P21912
GeneID	6390;
SDS-PAGE MW	31kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa, 293 cells using SDHB Polyclonal Antibody.

Note: This product is for in vitro research use only