

MCM6 Polyclonal Antibody

Catalog No: #41827

Package Size: #41827-1 50ul #41827-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MCM6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	MCM6 Polyclonal Antibody detects endogenous levels of MCM6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human MCM6 .
Target Name	MCM6
Other Names	MCM6; DNA replication licensing factor MCM6; p105MCM
Accession No.	Swiss-Prot: Q14566NCBI Gene ID: 4175
Uniprot	Q14566
GeneID	4175;
SDS-PAGE MW	90kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using MCM6 Polyclonal Antibody .

Note: This product is for in vitro research use only