

SNRPN Polyclonal Antibody

Catalog No: #41830

Package Size: #41830-1 50ul #41830-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	SNRPN Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	SNRPN Polyclonal Antibody detects endogenous levels of SNRPN protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human SNRPN.
Target Name	SNRPN
Other Names	SNRPN; HCERN3; SMN; Small nuclear ribonucleoprotein-associated protein N; snRNP-N; Sm protein D; Sm-D; Sm protein N; Sm-N; SmN; Tissue-specific-splicing protein
Accession No.	Swiss-Prot: P63162NCBI Gene ID: 6638
Uniprot	P63162
GeneID	6638;8926;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using SNRPN Polyclonal Antibody.

Note: This product is for in vitro research use only