

# Bovine Pyrroline-5-carboxylate reductase 2 (PYCR2) ELISA Kit



Catalog No: #EK7906

Orders: order@signalwayantibody.com

Package Size: #EK7906-1 48T #EK7906-2 96T

Support: tech@signalwayantibody.com

## Description

Product Name	Bovine Pyrroline-5-carboxylate reductase 2 (PYCR2) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Bovine (Bos taurus; Cattle)
Other Names	RP4-559A3.4; FLJ54750; P5CR2; pyrroline 5-carboxylate reductase isoform
Accession No.	Q17QJ7
Uniprot	Q17QJ7
GeneID	504987;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

## Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

**Detection Method:**SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate PYCR2 in samples. An antibody specific for PYCR2 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyPYCR2 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PYCR2 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PYCR2 bound in the initial step. The color development is stopped and the intensity of the color is measured.**Product Overview:**PYCR2, Belongs to the pyrroline-5-carboxylate reductase family.

Pyrroline-5-carboxylate reductase(P5CR) is the enzyme that catalyzes the terminal step in the biosynthesis of proline from glutamate, the NAD(P) dependent oxidation of 1-pyrroline-5-carboxylate into proline.

Pyrroline-5-carboxylate reductase, which converts pyrroline-5-carboxylate to proline, has been identified in human erythrocytes. The level of pyrroline-5-carboxylate reductase activity in these cells is comparable to the activity levels of major erythrocyte enzymes. The physiologic function of the enzyme in erythrocytes cannot be related to its function in other tissues, i.e., producing proline for protein synthesis.

---

Note: This product is for in vitro research use only