

IL-18R Polyclonal Antibody

Catalog No: #41844

Package Size: #41844-1 50ul #41844-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-18R Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IL-18R Polyclonal Antibody detects endogenous levels of IL-18R protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-18R.
Target Name	IL-18R
Other Names	IL18R1; IL1RRP; Interleukin-18 receptor 1; IL-18R-1; IL-18R1; CD218 antigen-like family member A; CDw218a; IL1 receptor-related protein; IL-1Rrp; IL1R-rp; CD218a
Accession No.	Swiss-Prot: Q13478NCBI Gene ID: 8809
Uniprot	Q13478
GeneID	8809;
SDS-PAGE MW	63kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of 293 cells using IL-18R Polyclonal Antibody.

Note: This product is for in vitro research use only