

NKp44 Polyclonal Antibody

Catalog No: #41852

Package Size: #41852-1 50ul #41852-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	NKp44 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	NKp44 Polyclonal Antibody detects endogenous levels of NKp44 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human NKp44.
Target Name	NKp44
Other Names	NCR2; LY95; Natural cytotoxicity triggering receptor 2; Lymphocyte antigen 95 homolog; NK cell-activating receptor; Natural killer cell p44-related protein; NK-p44; NKp44; CD336
Accession No.	Swiss-Prot: O95944NCBI Gene ID: 9436
Uniprot	O95944
GeneID	9436;
SDS-PAGE MW	44kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of PC-12 cells using NKp44 Polyclonal Antibody.

Note: This product is for in vitro research use only