

IL-33 Polyclonal Antibody

Catalog No: #41858

Package Size: #41858-1 50ul #41858-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-33 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC-p
Species Reactivity	Human,Mouse,Rat
Specificity	IL-33 Polyclonal Antibody detects endogenous levels of IL-33 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-33.
Target Name	IL-33
Other Names	IL33; C9orf26; IL1F11; NFHEV; Interleukin-33; IL-33; Interleukin-1 family member 11; IL-1F11; Nuclear factor from high endothelial venules; NF-HEV
Accession No.	Swiss-Prot: O95760NCBI Gene ID: 90865
Uniprot	O95760
GeneID	90865;
SDS-PAGE MW	31kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using IL-33 Polyclonal Antibody.

Note: This product is for in vitro research use only