

POLR3G Polyclonal Antibody

Catalog No: #41860

Package Size: #41860-1 50ul #41860-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	POLR3G Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms
Specificity	POLR3G Polyclonal Antibody detects endogenous levels of POLR3G protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human POLR3G.
Target Name	POLR3G
Other Names	POLR3G; DNA-directed RNA polymerase III subunit RPC7; RNA polymerase III subunit C7; DNA-directed RNA polymerase III subunit G; RNA polymerase III 32 kDa subunit; RPC32
Accession No.	Swiss-Prot: O15318NCBI Gene ID: 10622
Uniprot	O15318
GeneID	10622;
SDS-PAGE MW	26kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images



Western Blot analysis of PC12 cells using POLR3G Polyclonal Antibody.

Note: This product is for in vitro research use only