

## UBA1 Polyclonal Antibody

Catalog No: #41864

Package Size: #41864-1 50ul #41864-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	UBA1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	UBA1 Polyclonal Antibody detects endogenous levels of UBA1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human UBA1.
Target Name	UBA1
Other Names	UBA1; A1S9T; UBE1; Ubiquitin-like modifier-activating enzyme 1; Protein A1S9; Ubiquitin-activating enzyme E1
Accession No.	Swiss-Prot: P22314NCBI Gene ID: 7317
Uniprot	P22314
GeneID	7317;
SDS-PAGE MW	118kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of PC12 cells using UBA1 Polyclonal Antibody.

---

Note: This product is for in vitro research use only