

EF-1 α 1/2 Polyclonal Antibody

Catalog No: #41865

Package Size: #41865-1 50ul #41865-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	EF-1 α 1/2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	EF-1 α 1/2 Polyclonal Antibody detects endogenous levels of EF-1 α 1/2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human EF-1 α 1/2.
Target Name	EF-1 α 1/2
Other Names	EEF1A1; EEF1A; EF1A; LENG7; Elongation factor 1-alpha 1; EF-1-alpha-1; Elongation factor Tu; EF-Tu; Eukaryotic elongation factor 1 A-1; eEF1A-1; Leukocyte receptor cluster member 7; EEF1A2; EEF1AL; STN; Elongation factor 1-alpha 2; EF-1-alpha-2; Eukaryoti
Accession No.	Swiss-Prot: P68104/Q05639/Q5VTE0NCBI Gene ID: 1915/1917
Uniprot	P68104
GeneID	1915;
SDS-PAGE MW	50kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using EF-1 α 1/2 Polyclonal Antibody

Note: This product is for in vitro research use only