

Human Prolactin-releasing hormone (PRH) ELISA Kit

Catalog No: #EK8228



Package Size: #EK8228-1 48T #EK8228-2 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human Prolactin-releasing hormone (PRH) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human (Homo sapiens)
Other Names	PRH; PRRP; preprolactin-releasing peptide
Accession No.	P81277
Uniprot	P81277
GeneID	51052;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate PRLH in samples. An antibody specific for PRLH has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyPRLH present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PRLH is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PRLH bound in the initial step. The color development is stopped and the intensity of the color is measured.**Product Overview:**Prolactin-releasing hormone, also known as PRLH, is a hypothetical human hormone or hormone releasing factor. Existence of this factor has been hypothesized as prolactin is the only currently known hormone for which almost exclusively negative regulating factors are known (such as dopamine, leukemia inhibitory factor, some prostaglandins) but few stimulating factors. While many prolactin stimulating and enhancing factors are well known (such as thyrotropin-releasing hormone, oxytocin, vasoactive intestinal peptide and estrogen) those have primary functions other than stimulating prolactin release and the search for hypothetical releasing factor or factors continues.

The prolactin-releasing peptide identified in 1998 was a candidate for this function, however as of 2008 it appears its function is not yet completely elucidated.

Note: This product is for in vitro research use only