

CD47 Polyclonal Antibody

Catalog No: #41879

Package Size: #41879-1 50ul #41879-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	CD47 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF(P) IHC(P) ELISA
Species Reactivity	Hu Rt Ms
Specificity	CD47 Polyclonal Antibody detects endogenous levels of CD47 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human CD47.
Target Name	CD47
Other Names	CD47; MER6; Leukocyte surface antigen CD47; Antigenic surface determinant protein OA3; Integrin-associated protein; IAP; Protein MER6; CD47
Accession No.	Swiss-Prot: Q08722NCBI Gene ID: 961
Uniprot	Q08722
GeneID	961;
SDS-PAGE MW	35kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.


Not yet tested in other applications.

Images



Western Blot analysis of 3T3 cells using CD47 Polyclonal Antibody

Western Blot analysis of NIH-3T3 cells using CD47 Polyclonal Antibody



Note: This product is for in vitro research use only