

GMF- β Polyclonal Antibody

Catalog No: #41883

Package Size: #41883-1 50ul #41883-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	GMF- β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	GMF- β Polyclonal Antibody detects endogenous levels of GMF- β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human GMF- β .
Target Name	GMF- β
Other Names	GMFB; Glia maturation factor beta; GMF-beta
Accession No.	Swiss-Prot: P60983NCBI Gene ID: 2764
Uniprot	P60983
GeneID	2764;
SDS-PAGE MW	17kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

Images

Western Blot analysis of HeLa cells using GMF- β Polyclonal Antibody

Note: This product is for in vitro research use only