

Bovine Putative histone-lysine N-methyltransferase PRDM6 (PRDM6) ELISA Kit



Catalog No: #EK8347

Orders: order@signalwayantibody.com

Package Size: #EK8347-1 48T #EK8347-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Bovine Putative histone-lysine N-methyltransferase PRDM6 (PRDM6) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Bovine (Bos taurus; Cattle)
Other Names	PR-domain zinc finger protein 6
Accession No.	A6QPM3
Uniprot	A6QPM3
GeneID	519857;
Storage	<p>The stability of ELISA kit is determined by the loss rate of activity. The loss rate of this kit is less than 5% within the expiration date under appropriate storage condition.</p> <p>The loss rate was determined by accelerated thermal degradation test. Keep the kit at 37C for 4 and 7 days, and compare O.D.values of the kit kept at 37C with that of at recommended temperature. (referring from China Biological Products Standard, which was calculated by the Arrhenius equation. For ELISA kit, 4 days storage at 37C can be considered as 6 months at 2 - 8C, which means 7 days at 37C equaling 12 months at 2 - 8C).</p>

Application Details

Detect Range:Request Information

Sensitivity:Request Information

Sample Type:Serum, Plasma, Other biological fluids

Sample Volume: 1-200 µL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:SandwichTest principle:This assay employs a two-site sandwich ELISA to quantitate PRDM6 in samples. An antibody specific for PRDM6 has been pre-coated onto a microplate. Standards and samples are pipetted into the wells and anyPRDM6 present is bound by the immobilized antibody. After removing any unbound substances, a biotin-conjugated antibody specific for PRDM6 is added to the wells. After washing, Streptavidin conjugated Horseradish Peroxidase (HRP) is added to the wells. Following a wash to remove any unbound avidin-enzyme reagent, a substrate solution is added to the wells and color develops in proportion to the amount of PRDM6 bound in the initial step. The color development is stopped and the intensity of the color is measured.

Product Overview:PRDM6, is enriched in flk1(+) hematovascular precursor cells using a microarray-based approach. PRDM6 is also detected in embryonic and adult-derived endothelial cell lines. PRDM6 is co-localized with histone H4 and methylates H4-K20 (but not H3) in vitro and in vivo. Overexpression of PRDM6 in mouse embryonic endothelial cells induces apoptosis by activating caspase-3 and inducing G1 arrest. PRDM6 inhibits cell proliferation as determined by BrdU incorporation in endothelial cells, but not in rat aortic smooth muscle cells. PRDM6 is expressed by vascular precursors, has differential effects in endothelial cells and smooth muscle cells, and may play a role in vascular precursor differentiation and survival by modulating local chromatin-remodeling activity within hematovascular subpopulations during development.

Note: This product is for in vitro research use only