

INDOL1 Polyclonal Antibody

Catalog No: #41885

Package Size: #41885-1 50ul #41885-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	INDOL1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	INDOL1 Polyclonal Antibody detects endogenous levels of INDOL1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human INDOL1.
Target Name	INDOL1
Other Names	IDO2; INDOL1; Indoleamine 2,3-dioxygenase 2; IDO-2; Indoleamine 2,3-dioxygenase-like protein 1; Indoleamine-pyrrole 2,3-dioxygenase-like protein 1
Accession No.	Swiss-Prot: Q6ZQW0NCBI Gene ID: 169355
Uniprot	Q6ZQW0
GeneID	169355;
SDS-PAGE MW	45kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa HepG2 cells using INDOL1 Polyclonal Antibody

Note: This product is for in vitro research use only