

TRPC3 Polyclonal Antibody

Catalog No: #41890

Package Size: #41890-1 50ul #41890-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	TRPC3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	TRPC3 Polyclonal Antibody detects endogenous levels of TRPC3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human TRPC3.
Target Name	TRPC3
Other Names	TRPC3; TRP3; Short transient receptor potential channel 3; TrpC3; Transient receptor protein 3; TRP-3; hTrp-3; hTrp3
Accession No.	Swiss-Prot: Q13507NCBI Gene ID: 7222
Uniprot	Q13507
GeneID	7222;
SDS-PAGE MW	97kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using TRPC3 Polyclonal Antibody

Note: This product is for in vitro research use only