

Netrin-1 Polyclonal Antibody

Catalog No: #41892

Package Size: #41892-1 50ul #41892-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Netrin-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu Ms Rt
Specificity	Netrin-1 Polyclonal Antibody detects endogenous levels of Netrin-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Netrin-1.
Target Name	Netrin-1
Other Names	NTN1; NTN1L; Netrin-1; Epididymis tissue protein Li 131P
Accession No.	Swiss-Prot: O95631NCBI Gene ID: 9423
Uniprot	O95631
GeneID	9423;
SDS-PAGE MW	66kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Netrin-1 Polyclonal Antibody

Note: This product is for in vitro research use only