

Prolactin Polyclonal Antibody

Catalog No: #41894

Package Size: #41894-1 50ul #41894-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Prolactin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	Prolactin Polyclonal Antibody detects endogenous levels of Prolactin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Prolactin.
Target Name	Prolactin
Other Names	PRL; Prolactin; PRL
Accession No.	Swiss-Prot: P01236NCBI Gene ID: 5617
Uniprot	P01236
GeneID	5617;
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of SKOV3 cells using Prolactin Polyclonal Antibody

Note: This product is for in vitro research use only