

G-CSFR Polyclonal Antibody

Catalog No: #41896

Package Size: #41896-1 50ul #41896-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	G-CSFR Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	G-CSFR Polyclonal Antibody detects endogenous levels of G-CSFR protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human G-CSFR.
Target Name	G-CSFR
Other Names	CSF3R; GCSFR; Granulocyte colony-stimulating factor receptor; G-CSF receptor; G-CSF-R; CD114
Accession No.	Swiss-Prot: Q99062NCBI Gene ID: 1441
Uniprot	Q99062
GeneID	1441;
SDS-PAGE MW	92kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of SKOV3 cells using G-CSFR Polyclonal Antibody

Note: This product is for in vitro research use only