

IFITM2 Polyclonal Antibody

Catalog No: #41900

Package Size: #41900-1 50ul #41900-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	IFITM2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	IFITM2 Polyclonal Antibody detects endogenous levels of IFITM2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human IFITM2.
Target Name	IFITM2
Other Names	IFITM2; Interferon-induced transmembrane protein 2; Dispanin subfamily A member 2c; DSPA2c; Interferon-inducible protein 1-8D
Accession No.	Swiss-Prot: Q01629NCBI Gene ID: 10581
Uniprot	Q01629
GeneID	10581;
SDS-PAGE MW	15kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using IFITM2 Polyclonal Antibody

Note: This product is for in vitro research use only