

## Neurotensin Polyclonal Antibody

Catalog No: #41902

Package Size: #41902-1 50ul #41902-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	Neurotensin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Hu Ms Rt
Specificity	Neurotensin Polyclonal Antibody detects endogenous levels of Neurotensin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Neurotensin.
Target Name	Neurotensin
Other Names	NTS; Neurotensin/neuromedin N
Accession No.	Swiss-Prot: P30990NCBI Gene ID: 4922
Uniprot	P30990
GeneID	4922;
SDS-PAGE MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western Blot: 1/500 - 1/2000.

## Images



Western Blot analysis of HELA cells using Neurotensin Polyclonal Antibody

---

Note: This product is for in vitro research use only