

ANG I Polyclonal Antibody

Catalog No: #41909

Package Size: #41909-1 50ul #41909-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	ANG I Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Hu
Specificity	ANG I Polyclonal Antibody detects endogenous levels of ANG I protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human ANG I.
Target Name	ANG I
Other Names	ANG; RNASE5; Angiogenin; Ribonuclease 5; RNase 5
Accession No.	Swiss-Prot: P03950NCBI Gene ID: 283
Uniprot	P03950
GeneID	283;
SDS-PAGE MW	17kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of 293T Hela 823 cells using ANG I Polyclonal Antibody

Note: This product is for in vitro research use only